

Configuration & Paths

Table of contents

1 Overview	2
2 Data flow	2
3 Setup instructions	3
3.1 Step 1 — Clone the repository	3
3.2 Step 2 — Create the local data folder	3
3.3 Step 3 — Edit the config files	4
3.4 Step 4 — Add MATLAB paths	4
3.5 Step 5 — Install the Python environment	4
4 Config file reference	5
4.1 MATLAB — config/get_config.m	5
4.2 Python — config/config.py	6
5 Network drive structure	7
6 Local data folder structure	8
7 Automation scripts	8
8 Troubleshooting	9

Quick reference

Both MATLAB and Python have a single editable path (`project_root / PROJECT_ROOT`) that you set once per computer. All other paths are derived automatically.

1 Overview

The freely-walking-optomotor system runs across three computers, each with a specific role. A centralised configuration in the code repository ensures all scripts find data in the correct locations.

Computer	Role	Config fields used
Acquisition rig (Windows)	Runs protocols, records video	<code>rig_data_folder</code> , <code>bias_config</code> , <code>SOURCE_ROOT</code>
Processing machine (Windows)	Automated tracking & processing	<code>project_root</code> + all local paths + all <code>NETWORK_*</code> paths
Analysis computer (Mac/any)	Manual analysis, plotting, dashboard	<code>project_root</code> + local paths only

The two config files live in the `config/` directory at the repository root:

- **MATLAB:** `config/get_config.m` — returns a struct via `cfg = get_config()`
- **Python:** `config/config.py` — exports constants via `from config.config import ...`

2 Data flow

The automation pipeline moves data through four stages:



1. **Acquisition rig** records video and LOG files to `rig_data_folder`
2. `monitor_and_copy` transfers completed folders to the network drive (`0_unprocessed/`)
3. `monitor_and_track` copies to local, runs FlyTracker, moves tracked data to `1_tracked/`
4. `daily_processing` runs MATLAB processing, saves results and figures, moves data to `2_processed/`

The processing machine maintains a local copy of the data under `project_root` and synchronises with the network drive. Analysis computers can work from either a local copy or a mounted network drive.

3 Setup instructions

3.1 Step 1 — Clone the repository

```
git clone https://github.com/leburnett/freely-walking-optomotor.git
```

3.2 Step 2 — Create the local data folder

Create the following folder structure (or point to an existing one):

```
your_data_root/  
  DATA/  
    00_unprocessed/  
    01_tracked/  
    02_processed/
```

```
results/  
figures/
```

3.3 Step 3 — Edit the config files

MATLAB — edit `config/get_config.m` (one line):

```
cfg.project_root = '/path/to/your_data_root';
```

Python — edit `config/config.py` (one line):

```
PROJECT_ROOT = Path("/path/to/your_data_root")
```

Both files contain commented examples for Windows and Mac.

3.4 Step 4 — Add MATLAB paths

Run `setup_path.m` once per MATLAB session (or add it to your `startup.m`):

```
run('/path/to/freely-walking-optomotor/setup_path.m')
```

This adds all `src/` subdirectories to the MATLAB path so that functions like `get_config()`, `process_freely_walking_data()`, and the protocol helper functions are available.

3.5 Step 5 — Install the Python environment

```
cd freely-walking-optomotor/python/freely-walking-python  
pixi install
```

This installs all Python dependencies (numpy, pandas, scipy, dash, etc.) into an isolated environment managed by [pixi](#).

4 Config file reference

4.1 MATLAB — config/get_config.m

Call `cfg = get_config()` to get the configuration struct.

Editable:

Field	Description
<code>cfg.project_root</code>	Root of your local data folder — edit this per computer

Derived from `project_root` (processing machine + analysis computer):

Field	Path	Used by
<code>cfg.data_unprocessed</code>	DATA/00_unprocessed/	Tracking pipeline
<code>cfg.data_tracked</code>	DATA/01_tracked/	<code>process_freely_walking_data</code>
<code>cfg.data_processed</code>	DATA/02_processed/	Archive
<code>cfg.results</code>	results/	Processing, analysis, plotting, dashboard
<code>cfg.figures</code>	figures/	Overview figure output

Derived from `repo_root` (all machines):

Field	Path	Used by
<code>cfg.repo_root</code>	Auto-detected git repo root	Internal
<code>cfg.patterns</code>	src/patterns/Patterns_protocols/	Protocols, pattern tools
<code>cfg.calibration_file</code>	src/tracking/calibration_data/	<code>batch_track_umf</code>

Acquisition rig only:

Field	Value	Used by
<code>cfg.rig_data_folder</code>	C:\MatlabRoot\FreeWalkOptimization\data	<code>optical_video_and_folders</code> (protocols)
<code>cfg.bias_config</code>	C:\MatlabRoot\...\bias_config\free_walk	<code>optical_video_and_folders</code> (BIAS camera)

Network drive:

Field	Value	Used by
cfg.group_drive	smb://prfs.hhmi.org/reislab/data/cokey/data/	Access from Mac

4.2 Python — config/config.py

Import variables with `from config.config import DATA_TRACKED, RESULTS_PATH`.

Editable:

Variable	Description
PROJECT_ROOT	Root of your local data folder — edit this per computer

Derived from PROJECT_ROOT (processing machine + analysis computer):

Variable	Path	Used by
DATA_UNPROCESSED	DATA/00_unprocessed/	monitor_and_track
DATA_TRACKED	DATA/01_tracked/	daily_processing
DATA_PROCESSED	DATA/02_processed/	daily_processing
RESULTS_PATH	results/	daily_processing, dashboard
FIGURES_PATH	figures/	daily_processing

Derived from REPO_ROOT:

Variable	Path	Used by
REPO_ROOT	Auto-detected repo root	Automation scripts
PATTERNS_DIR	src/patterns/PatternsDocumentation/	Documentation generator
PROTOCOLS_DIR	src/protocols/	Documentation generator

Network drive paths (processing machine only):

Variable	Path	Used by
NETWORK_ROOT	\\prfs.hhmi.org\reislab\data/cokey	Access to network paths

Variable	Path	Used by
NETWORK_UNPROCESSED	...\data\0_unprocessed	monitor_and_copy, monitor_and_track
NETWORK_TRACKED	...\data\1_tracked	monitor_and_track, daily_processing
NETWORK_PROCESSED	...\data\2_processed	daily_processing
NETWORK_RESULTS	...\exp_results	daily_processing
NETWORK_FIGS	...\exp_figures\overview_figs	daily_processing

Acquisition rig (rig computer only):

Variable	Value	Used by
SOURCE_ROOT	C:\MatlabRoot\FreeWalk	monitor_and_copy

5 Network drive structure

The group network drive is the central archive for all experiment data:

```

\\prfs.hhmi.org\reiserlab\oaky-cokey\
data\
  0_unprocessed/  <- raw data from rig (via monitor_and_copy)
  1_tracked/     <- tracked data (via monitor_and_track)
  2_processed/   <- processed data with MP4 videos (via daily_processing)
exp_results/     <- .mat result files (via daily_processing)
exp_figures\
  overview_figs/ <- overview PDF/PNG figures (via daily_processing)

```

i UNC vs SMB paths

The Python config uses Windows UNC format (\\server\share) while the MATLAB config uses SMB format (smb://server/share/data/). The Python NETWORK_ROOT points to the share root (oaky-cokey); the MATLAB `cfg.group_drive` points to the `data/` subdirectory underneath.

6 Local data folder structure

Both the processing machine and analysis computers use the same folder layout under `project_root / PROJECT_ROOT`:

```
project_root/  
  DATA/  
    00_unprocessed/  <- raw data staged for tracking  
    01_tracked/      <- tracked data awaiting processing  
      YYYY_MM_DD/  
        protocol_name/  
          strain/  
            sex/  
              HH_MM_SS/  
                LOG_YYYYMMDD_HHMMSS.mat  
                REC_.ufmf  
                tracking/  
                  trx.mat  
                  *-feat.mat  
    02_processed/    <- fully processed archive  
  results/  
    protocol_27/  
      strain/  
        sex/  
          YYYYMMDD_HHMMSS_strain_protocol_sex_data.mat  
  figures/  
    overview_figs/  
    ...
```

i Local vs network numbering

Local folders use zero-padded numbers (00_, 01_, 02_) while the network drive uses single digits (0_, 1_, 2_). This is intentional — it makes it easy to tell whether you are looking at a local or network path.

7 Automation scripts

The three automation scripts in `python/automation/` run on the processing machine and move data through the pipeline:

Script	Runs on	Config imports	Purpose
<code>monitor_and_copy</code>	Acquisition rig	<code>SOURCE_ROOT</code> , <code>NETWORK_UNPROCESSED</code>	Watches for completed recordings, copies to network
<code>monitor_and_track</code>	Processing machine	<code>DATA_UNPROCESSED</code> , <code>DATA_TRACKED</code> , <code>NETWORK_*</code> , <code>REPO_ROOT</code>	Copies from network, runs FlyTracker, archives tracked data
<code>daily_processing</code>	Processing machine	<code>DATA_TRACKED</code> , <code>DATA_PROCESSED</code> , <code>NETWORK_*</code> , <code>RESULTS_PATH</code> , <code>FIGURES_PATH</code> , <code>REPO_ROOT</code>	Runs MATLAB processing, copies results to network

See the [Data Organisation](#) page for details on what each processing stage produces.

8 Troubleshooting

“Path not found” errors in MATLAB Run `setup_path.m` to add all `src/` subdirectories to the MATLAB path, then check that `get_config().project_root` points to an existing directory:

```
cfg = get_config();
assert(isfolder(cfg.project_root), 'project_root does not exist: %s', cfg.project_root);
```

Network paths unreachable The `NETWORK_*` paths require the network drive to be mounted. On Windows this is typically automatic via domain login. On Mac, connect via Finder: Go > Connect to Server > `smb://prfs.hhmi.org/reiserlab/oaky-cokey`.

Python import errors Make sure `sys.path.insert` reaches the repo root before importing from `config.config`. The depth depends on the script’s location in the directory tree. For example, scripts in `python/automation/daily_processing/` need four levels:

```
sys.path.insert(0, str(Path(__file__).parent.parent.parent.parent))
```

Dashboard can’t find data The dashboard reads from `RESULTS_PATH / "protocol_27"` by default. Ensure you have run `process_freely_walking_data` for the relevant dates, or point the dashboard to a different protocol directory.