

Data

Table of contents

1 Where to locate the data 1

[Link](#) to documentation on how the automatic logging system was set up.

[Link](#) to the google sheet with the automatic logging of experiments.

Each protocol run generates:

- **Video:** UFMF format via SimpleBiasCameraInterface
- **LOG.mat:** Trial timing, conditions, temperature, metadata
- **Google Sheets:** Optional export for experiment tracking

Output folder structure:

```
project_data_folder/  
  YYYY-MM-DD/  
    protocol_name/  
      strain/  
        sex/  
          HH-MM-SS/  
            LOG_YYYYMMDD_HHMMSS.mat  
            video.ufmf
```

1 Where to locate the data