

protocol_23

Table of contents

1	Protocol Parameters	1
2	Conditions	1
3	Pattern Previews	2
3.0.1	Pattern 45	2
3.0.2	Pattern 46	2

Protocol_23.m - test different contrast values for ON and OFF bars Based on protocol 19
Includes: - ALL OFF interval frame. - 8 pixel ON and OFF bar fixation - 180 deg apart —
varying contrast values (2, 4, and 6)

1 Protocol Parameters

Parameter	Value
Conditions	2
Patterns Used	2
t_acclim	20s
t_pause	0.015

2 Conditions

#	Pattern	Speed	Duration	Description
1	Pattern 45	9	-s	bar fixation - 16px ON
2	Pattern 46	9	-s	bar fixation - 16px OFF

3 Pattern Previews

3.0.1 [Pattern 45](#)



3.0.2 [Pattern 46](#)

